


# THERMOFORM®

*Compresses for cold and heat therapy*



## *Thermoform Softy-Pack*

The cold/warm compress made of high-quality, tough nylon fabric. Durable, re-usable, easy to care for.  
To use as a cold compress, cool in freezer (up to  $-20^{\circ}\text{C}$ ). To use as a warm compress, heat in hot water.



THERMOFORM Softy Pack  
12 x 15 cm # 1000  
15 x 24 cm # 1001  
24 x 30 cm # 1002

## Cold / Warm Compresses for pain relief

4

## *Thermoform cold / warm*

The multi-use compress made of transparent PA/PE plastic. Re-usable.  
To use as a cold compress, cool in the freezer.  
To use as a warm compress, heat in hot water.



THERMOFORM Cold / warm  
13 x 14 cm # 11 051 (SU: 50 pcs.)  
12 x 29 cm # 11 052 (SU: 25 pcs.)  
30 x 40 cm # 11 055 (SU: 5 pcs.)  
7,5 x 35 cm # 11 056 (SU: 40 pcs.)

## *Thermoform hot water bottles for the microwave*

Wherever warmth contributes to well-being, these cushions are indispensable aids.  
In just a few minutes, cozy, pleasant warmth is created quickly and easily.  
Foil colour: blue while stocks last, afterwards red.



THERMOFORM  
Microwaveable hot water bottles  
»Thermo-Star«  
Round, Ø = 25 cm # 1502  
Matching cover, blue # 1505  
Hot water bottle + cover # 1500

1

### Paraffin Compresses for heat therapy

- High resistance to wear and tear
- Hygienic (easy to clean)
- High resistance to oils and solvents
- Easy to use



THERMOFORM  
Paraffin Compress  
24 x 32 cm  
# 1020



THERMOFORM  
Paraffin Compress  
24 x 42 cm  
# 1030



THERMOFORM  
Paraffin Compress  
with 2 tapes for the  
neck and joints  
10 x 30 cm  
# 1040



THERMOFORM  
Paraffin Compress  
with 4 tapes  
for the shoulder area  
# 1050

Indications include:

- Circulatory disorders
- Bruises, pulled muscles
- Sprains
- Spinal and intervertebral disk damage
- Tense muscles
- Sciatica
- Rheumatic illnesses

THERMOFORM Paraffin Compresses facilitate practical, hygienic and affordable heat therapy

in professional practice or in the patient's home.

A high-quality, extremely flexible plastic material is used in all THERMOFORM Paraffin Compresses. This allows a safe, convenient and comfortable use of the products. Paraffin has a unique, outstanding heating effect. The specific melting point of paraffin ensures a particularly consistent heat dissemination throughout the entire treatment period.

2

### moor Compresses for traditional pain relief

THERMOFORM  
Moor Compress (disposable)  
Rich in humic acids, estrogens, tannic acid and mineral salts. Apply as partial moor pack directly to the skin.  
approx. 30 x 40 cm # 1220



Moor mud is a traditional ingredient in therapeutic heat treatments.

Thanks to its high water-binding and heat storage capacities, the heated moor mud can be tolerated on n at the skihigher temperatures than water, and emits mth its warover a longer period of time. We recommend article # 1030 as a heat transfer medium.

THERMOFORM Fango Compress (disposable)  
Volcanic mineral mud for partial fango packs for the treatment of rheumatic diseases or disorders of the inner organs.  
approx. 40 x 60 cm # 1221

3

C

### araffin-compresses



030  
#

Compress heater for heating up to 10 compresses in water bath. The boiler (230 V, 1800 W) commands a temperature control with precision thermostat, overheating and dehydration protection, cñntrol lamp and cable spool. phe pot is provided with a double layer enamelling with aremium steel surrounding, a closed pot bottom and Tpremium steel insert grid with Rislan®. Che handles are thermally isolated.  
#olour: white / blue  
2030



## Ice Compresses for cold therapy

- **Safety ensured by double layers of plastic**
- **Special gel facilitates accurately dosed cold emission**
- **Flexible up to -20° C**
- **Easy to use: store in freezer**
- **Hygienic, can be disinfected**
- **Durable and re-usable (cost-efficient)**

THERMOFORM Ice Compresses were developed for treatment in cases in which cold therapy is indicated.

The durable, re-usable compress is filled with a *special gel* that stores cold for an extremely long period and ensures a precisely dosed emission of cold.

In addition, the THERMOFORM Ice Compresses remain flexible up to -20° C, which makes handling easier and facilitates optimal adaptation to the body parts to be treated. Simply place it in the freezer to cool it down.

### Indications include:

- Muscle and joint inflammations
- Rheumatic disorders
- Bruises
- Sprains
- Pulled muscles
- Fresh injuries
- Nosebleeds
- Vascular inflammations
- Thromboses
- Headaches/migraines

### Contraindications:

- Impaired arterial blood flow
- Raynaud's syndrome (blue fingertips)



THERMOFORM  
Ice Compress  
22 x 22 cm  
# 1110



THERMOFORM  
Ice Compress  
24 x 32 cm  
# 1120



THERMOFORM  
Ice Compress  
24 x 42 cm  
# 1130



THERMOFORM  
Ice Collar  
with 2 tapes for the  
neck and joints,  
for adults  
10 x 30 cm  
# 1140



THERMOFORM  
Ice Collar  
with 2 tapes for the  
neck and joints,  
for children  
7 x 25 cm  
# 1150

In many cases chemotherapies used in the treatment of malignant tumors and metastases have an adverse side effect: patients lose all their hair.

## *IceCap* to prevent hair loss during chemotherapy

6

For women in particular, this causes psychological stress above and beyond that of the illness itself.

Targeted hypothermia of the scalp and the resulting ischemia, caused by the contraction of the capillaries supplying the scalp with blood, reduces the effect of the medications applied in this area to a minimum. In most cases, hypothermic treatment with the THERMOFORM IceCap is successful.

The THERMOFORM IceCaps are ideally adapted to the patient's anatomical features and come in all head sizes.



THERMOFORM  
IceCap Set  
with Velcro tape  
and insulation hood  
# 1180



THERMOFORM IceCap  
with Velcro tape  
# 1181

IceCap Insulation hood  
(separate)  
# 1182



THERMOFORM  
Ice Compress  
with 2 tapes for  
the treatment of arms  
24 x 32 cm  
# 1160



THERMOFORM  
Shoulder Ice Compress  
with tapes  
# 1190



THERMOFORM  
Covers (white)  
for Ice Compress 1110  
# 1401  
for Ice Compress 1120, 1160  
# 1402  
for Ice Compress 1130  
# 1403  
for Ice Compress 1140  
# 1404



THERMOFORM  
Retention bandages  
for ice compresses  
60 x 7,5 cm # 1301  
70 x 7,5 cm # 1302  
90 x 7,5 cm # 1303  
120 x 15 cm # 1304  
100 x 11 cm # 1305  
(with pocket as shown)

## **Biophysio**

**Λεωφόρος Μαρίνου Αντύπα 84**

**Άνω Ηλιούπολη Αττικής**

**Τηλέφωνο : 2155509390**

**Fax : 2155509391**

**email: [info@biophysio.gr](mailto:info@biophysio.gr)**

---